
# Mathematica code for AFM-IR Ratiomap (v2.1)

Contact: jongcheol1422@gmail.com



### v2.1 updates

1. Defined a function that generate the 'ratio' table. You don't need to copy&paste a few lines that divide a table by another table manuall 1-1. How to use? You need to run & save the results to a table with a

proper name.

1164/1745 --> [S2-S1], great @ ratio[rIR1164b, rIR1745] ir1164ov1745 = temp; IRmap[ir1164ov1745, 20] ir1164ov1745 Intensity 0.5 1. 1.5 2. 2.5 3. 3.5 4. 4.5 5. 70.2 62.4

2. Defined a new function for the Contour Plot. Now, you can freely control the color scheme range using 'IRmap2' function.

2-2. How to use? 4 arguments. [table, min, max, bin]

```
IRmap2[ir1164ov1745, 1.9, 7, 20] (*wvn,min,max,bin*)
                          ir1164ov1745
                                                                   Intensity
                                                                      6.49
                                                                      5.98
   4.5
                                                                      5.47
                                                                      4.96
                                                                      4.45
                                                                      3.94
                                                                      -3.43
   3.5
                                                                      2.92
                                                                      2.41
```

**b.** Normalization - Ratio map **>** 

cut1 = 0; (\*cut out outliers\*)

```
cut2 = 1000; (*cut out outliers*)
ratio[aa, bb]:= Module[{},
  temp = Table[0, {i, 1, dim}, {j, 1, dim}];
  For [i = 1, i \le dim,
   For [j = 1, j \le dim,
    If[cut1 < aa[i, j] / bb[i, j] < cut2,</pre>
     temp[i, j] = aa[i, j] / bb[i, j];
    j++];
                                c. Imaging D
   i++]
                                 0. inputs 3
```

```
(*Unnormalized*)
SetAttributes[IRmap, HoldFirst];
IRmap[temp_, cont_] := Module[{name1},
 name1 = SymbolName[Unevaluated[temp]];
 ListContourPlot[Evaluate[temp], PlotRange → All, Frame
   AspectRatio → ny / nx, ColorFunction → "Rainbow", Conto
   PlotLegends → Placed[BarLegend[Automatic, LegendMarke
```

```
(*Normalized*)
ClearAll[IRmap2, min, max, dd];
SetAttributes[IRmap2, HoldFirst];
IRmap2[temp_, min_, max_, dd_] := Module[{name2, leg},
  name2 = SymbolName[Unevaluated[temp]];
  leg = Round[Table[i, {i, min, max, (max - min) / dd}], 0.0
  ListContourPlot [Evaluate[temp], PlotRange → All, Frame
   AspectRatio \rightarrow \frac{ny}{nx}, GridLines \rightarrow {Range[0, nx, dim/Lx * o
   ColorFunction → Function[{z}, ColorData["Rainbow"][Re
   ColorFunctionScaling → False,
   Contours → leg, ContourStyle → None, (* Related to col
   PlotLegends → Placed[BarLegend[Automatic, LegendMa
   PlotLabel → Dynamic[name2]
(*PlotLegends->Placed[BarLegend[{"Rainbow",{min,max}}
```

### v2.1 updates

3. Make sure that the imported data files are properly assigned in the later step. To do that, Print function was added.

#### 1. Data import

```
allFiles = FileNames["*.txt", "sample data"]; (*import .txt files in the directory name 'data'*)
allData = Import[#, "Data"] & /@ allFiles;
allData // Dimensions(*Check how many .txt files were imported*)
Print["Imported file names: ", allFiles](* Make sure to chekc the order of the file imported! *)
```

{6, 65537, 1}


Imported file names: {sample data/1034.txt, sample data/1060.txt, sample data/1151.txt, sample data/r1034.txt, sample data/r1060.txt

#### 3. Data name assignment


```
(*name each dataset for convenience*)
IR1034 = e[1];
IR1060 = e[2];
IR1151 = e[3];
rIR1034 = e[4];
rIR1060 = e[5];
rIR1151 = e[6];
```

# Mathematica code for AFM-IR Ratiomap (v2)


Contact: jongcheol1422@gmail.com



How to expoert Data from Bruker - Nanoscope software.



## Check if the txt files are below. First line is the value name. spl x spl values are in one column. Ver. 2



## Cells with light yellow background $\rightarrow$ You need to input proper numbers about your data information.

Explanations are in (\*exp\*)

Ver.2 (nx x ny data in one column, in .txt)

```
In[62]:= $Version
     ClearAll["Global`*"]
     SetDirectory[NotebookDirectory[]]; (*Mathematica file in the proper location*)
Out[62]= 12.3.0 for Mac OS X x86 (64-bit) (May 10, 2021)
   a. Data import and process D
     0. Data information
In[65]:= Lx = 5; (*um-scanned area *)
       Ly = 5; (*um-scanned area *)
       dim = 256; (*sample per line *)
       nx = dim;
       ny = dim;
     1. Data import
       allFiles = FileNames["*.txt", "sample data"]; (*import .txt files in the directory name 'data'*)
       allData = Import[#, "Data"] & /@ allFiles;
       allData // Dimensions (*Check how many .txt files were imported*)
Out[71]= \{10, 65537, 1\}
       allData2 = Table[i, {i, 1, Length[allData]}];
       For[i = 1, i ≤ Length[allData],
       allData2[i] = Drop[allData[i], 1]; (*Drop the first line (names)*)
       i++]
       allData // Dimensions
```