




# Mathematica code for AFM-IR Ratiomap (v2)

Contact: [jongcheol1422@gmail.com](mailto:jongcheol1422@gmail.com)




# How to export Data from Bruker – Nanoscope software.




Native → raw data (for Amp2)

Display → numbers after process in Nanoscope (for height)

LSB → R



**Check if the txt files are below. First line is the value name. spl x spl values are in one column. Ver. 2**



# Cells with light yellow background → You need to input proper numbers about your data information.

Explanations are in (\*exp\*)

Ver.2 (nx x ny data in one column, in .txt)

```
In[62]:= $Version
ClearAll["Global`*"]
SetDirectory[NotebookDirectory[]]; (*Mathematica file in the proper location*)

Out[62]= 12.3.0 for Mac OS X x86 (64-bit) (May 10, 2021)
```

## a. Data import and process

### 0. Data information

```
In[65]:= Lx = 5; (*um-scanned area *)
Ly = 5; (*um-scanned area *)
dim = 256; (*sample per line *)
nx = dim;
ny = dim;
```

### 1. Data import

```
In[70]:= allFiles = FileNames["*.txt", "sample data"]; (*import .txt files in the directory name 'data'*)
allData = Import[#, "Data"] & /@ allFiles;
allData // Dimensions (*Check how many .txt files were imported*)
```

```
Out[71]= {10, 65537, 1}
```

```
In[72]:= allData2 = Table[i, {i, 1, Length[allData]}];
For[i = 1, i ≤ Length[allData],
  allData2[[i]] = Drop[allData[[i]], 1]; (*Drop the first line (names)*)
  i++]
allData // Dimensions
```